

# uPA antibody

Catalog Number: 116622

## Product name

uPA antibody

WB: 1:200-1:2000

IP: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:20-1:200

## Antibody description

uPA Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human skin tissue, human lung tissue, human ovary tissue. Positive WB detected in HeLa cells, A375 cells, human brain tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 54 kDa

## Validations



## Preparation

This antibody was obtained by immunization of uPA recombinant protein (Accession Number: XM\_047425353). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116622(uPA antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IP, IF

## Dilutions



Recommended Dilution:



Immunohistochemical of paraffin-embedded human skin using Catalog No:116622(uPA antibody) at dilution of 1:50 (under 40x lens)

# uPA antibody

Catalog Number: 116622



Immunofluorescent analysis of HeLa cells using Catalog No:116622(uPA Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

IP Result of anti-uPA (IP:Catalog No:116622, 4ug; Detection:Catalog No:116622 1:1000) with HeLa cells lysate 2800ug.