

# TH antibody



Catalog Number: 107622

## Product name

TH antibody

## Specificity

Human; other species not tested.

## Antibody description

TH Mouse Monoclonal antibody. Positive WB detected in SH-SY5Y cells, PC-12 cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of TH recombinant protein (Accession Number: NM\_000360). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107622(TH Antibody) at dilution of 1:2000