

CCL4 Polyclonal Antibody



Catalog Number: 163322

Product name

CCL4 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CCL4

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 30 to the C-terminus of human CCL4 (NP_002975.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - CCL4 Polyclonal Antibody

Western blot analysis of extracts of SGC-7901 cells, using CCL4 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - CCL4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using CCL4 antibody at dilution of 1:100 (40x lens).