

Anti-HDAC2 (2D9) Mouse antibody



Catalog Number: 168322

Product name

Anti-HDAC2 (2D9) Mouse antibody

Ig Type

Mouse IgG2b

Specificity

Human, Mouse, Rat, Monkey

Applications

WB, ICC/IF

Antibody description

HDAC2 (2D9) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human HDAC2 protein fragments expressed in E.coli.

Validations

Formulation


PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of HDAC2 in THP-1, SHSY-5Y, K562, COS7, C6, 3T3 and HeLa cell lysates using HDAC2 mouse mAb(dilution 1:1000). Predicted band size: 60kDa. Observed band size: 60kDa.