

Anti-Cullin 5/CUL-5 Rabbit antibody



Catalog Number: 166722

Product name

Anti-Cullin 5/CUL-5 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

Cullin 5/CUL-5 Rabbit polyclonal antibody

Dilutions

WB: 1/1000-1/5000


Preparation

Antigen: Recombinant protein of human Cullin 5

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Western blot detection of Cullin 5/CUL-5 in C6,3T3,HeLa cell lysates using Cullin 5/CUL-5 Rabbit pAb(1:1000 diluted).Predicted band size:91kDa.Observed band size:91kDa.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal