

Anti-SFRP5 antibody



Catalog Number: 176722

Product name

Mouse IgG1

Anti-SFRP5 antibody

Applications

Specificity

WB, IP, IF

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1:100-1:1,000

Mouse monoclonal antibody to SFRP5

IP: 1-2 µg per 100-500 µg of total protein(1 ml of cell lysate)

Preparation

IF: 1:50-1:500

This antigen of this antibody was peptide

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 40%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C

Clonality

Monoclonal

Ig Type

Validations




Fig1: Western blot analysis of SARP3 expression in non-transfected (A) and human SARP3 transfected (B) 293T whole cell lysates.