

Anti-Caspase-1 Rabbit antibody



Catalog Number: 168722

Product name

Rabbit IgG

Anti-Caspase-1 Rabbit antibody

Applications

Specificity

WB, IP

Human

Dilutions

Antibody description

WB: 1/1000

Caspase-1 Rabbit monoclonal antibody

IP: 1/20


Preparation

Antigen: Recombinant protein of human Caspase-1

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Caspase-1 in HeLa cell lysates using Caspase-1 Rabbit mAb(1:500 diluted).Predicted band size:45kDa.Observed band size:45,20kDa.

Clonality

Monoclonal

Ig Type