

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 345) of the influenza hemagglutinin (A/reassortant/IVR-155(Victoria/210/2009 x Puerto Rico/8/1934)(H3N2))(ADI52838.1), termed as HA1, was expressed with a C-terminal polyhistidine tag.

Organism

H3N2

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The secreted recombinant HA1 subunit of influenza A H3N2

(A/reassortant/IVR-155(Victoria/210/2009 x Puerto Rico/8/1934)(H3N2)) comprises 340 amino acids and has a predicted molecular mass of 38.1 kDa. It migrates as an approximately 58-68 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N2 (A/reassortant/IVR-155) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE