

Anti-UGP2 antibody



Catalog Number: 175032

Product name

Anti-UGP2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Mouse monoclonal antibody to UGP2

Preparation

This antigen of this antibody was recombinant protein.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, FC

Dilutions

WB: 1:500

IHC-P: 1:50-1:100

FC: 1:50-1:100

Validations





Fig1: Western blot analysis of UGP2 on human liver tissue lysate using anti-UGP2 antibody at



1/200 dilution.

Fig2: Immunohistochemical analysis of paraffin-embedded rat skeletal muscle tissue using anti-UGP2 antibody. Counter stained with hematoxylin.




Fig3: Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-UGP2 antibody. Counter stained with hematoxylin.

Anti-UGP2 antibody

Catalog Number: 175032




Fig4: Immunohistochemical analysis of paraffin-embedded mouse heart tissue using anti-UGP2 antibody. Counter stained with hematoxylin.




Fig5: Flow cytometric analysis of LOVO cells with UGP2 antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti-mouse IgG was used as the secondary antibody.