

Anti-IL17RC antibody

Catalog Number: 176032

Product name

Anti-IL17RC antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to IL17RC

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human il17rc. aa range:721-770

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using IL-17RC Polyclonal Antibody