

# ZNF763 antibody



Catalog Number: 117032

---

## Product name

ZNF763 antibody

## Specificity

human; other species not tested.

## Antibody description

ZNF763 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF7 cells. Observed molecular weight by Western-blot: 46-47kd

## Preparation

This antibody was obtained by immunization of ZNF763 recombinant protein (Accession Number: NM\_001367172). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117032(ZNF763 antibody) at dilution of 1:1000