

Human LEFTY2/EBAF/LEFTA/LEFTYA (His tag) recombinant protein



Catalog Number: 519032

General Information

Gene Name Synonym

Left-right determination factor 2; Endometrial bleeding-associated factor; Left-right determination factor A; Protein lefty-2; Protein lefty-A; Transforming growth factor beta-4; TGF-beta-4; LEFTY2; EBAF; LEFTA; LEFTYA; TGFB4

Protein Construction

A DNA sequence encoding the Human LEFTY2/EBAF/LEFTA/LEFTYA Phe78-Pro366 is expressed fused with a Mouse Lefty-1 Propeptide (Leu22-Arg77) & 6His tag at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of PBS, pH7.4.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human LEFTY2/EBAF/LEFTA/LEFTYA (His tag)
recombinant protein