

RAB33A antibody

Catalog Number: 114432

Product name

RAB33A antibody

Specificity

human, mouse; other species not tested.

Antibody description

RAB33A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue. Positive IHC detected in human lymphoma. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of RAB33A recombinant protein (Accession Number: NM_004794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114432(RAB33A antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114432(RAB33A antibody) at dilution of 1:50 (under 10x lens)