

Human ACYP1/ACYPE (His tag) recombinant protein



Catalog Number: 518432

General Information

Gene Name Synonym

Acylphosphatase-1, ACYP1, Acylphosphatase, erythrocyte isozyme, Acylphosphatase, organ-common type isozyme, Acylphosphate phosphohydrolase 1, ACYPE

Protein Construction

A DNA sequence encoding the Human ACYP1/ACYPE Met1-Lys99 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 100mM NaCl, 20% Glycerol, 1mM DTT, pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ACYP1/ACYPE (His tag) recombinant protein