

# Galc antibody

Catalog Number: 110832

## Product name

Galc antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Galc Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A375 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 80kd, 30kd

## Preparation

This antibody was obtained by immunization of Galc recombinant protein (Accession Number: BC086671). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110832(Galc antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110832(Galc antibody) at dilution of 1:50 (under 10x lens)