

AQP5 Polyclonal Antibody



Catalog Number: 161832

Product name

AQP5 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to AQP5

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human AQP5 (NP_001642.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:200

Validations



Immunohistochemistry - AQP5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using AQP5 antibody at dilution of 1:50 (40x lens).