

# Human PPM1A/PPPM1A (His tag) recombinant protein



Catalog Number: 519832

## General Information

### Gene Name Synonym

Protein Phosphatase 1A, Protein Phosphatase 2C Isoform Alpha, PP2C-Alpha, Protein Phosphatase 1A, PPM1A, PPPM1A

### Protein Construction

A DNA sequence encoding the Human PPM1A/PPPM1A Gly2-Trp382 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PPM1A/PPPM1A (His tag) recombinant protein