

Dynamitin antibody

Catalog Number: 110132

Product name

Dynamitin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Dynamitin Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, human testis tissue, mouse lung tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_001261413). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110132(dynactin-2 antibody) at dilution of



1:800

IP Result of anti-dynactin-2 (IP:Catalog No:110132, 3ug; Detection:Catalog No:110132 1:800) with mouse brain tissue lysate 4000ug.