

# GPR58/TAAR2 antibody



Catalog Number: 111132

---

## Product name

GPR58/TAAR2 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

GPR58/TAAR2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse brain tissue, mouse cerebellum tissue. Observed molecular weight by Western-blot: 42 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of GPR58/TAAR2 recombinant protein (Accession Number: NM\_014626). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111132(TAAR2 Antibody) at dilution of 1:600

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality