

ANO1 Polyclonal Antibody



Catalog Number: 160532

Product name

ANO1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ANO1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 628-731 of human ANO1 (NP_060513.5).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:1000 - 1:2000


IF 1:100 - 1:200

Validations



Western blot - ANO1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ANO1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunofluorescence - ANO1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using ANO1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.