

AIRE Polyclonal Antibody



Catalog Number: 163532

Product name

AIRE Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to AIRE

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 30-130 of human AIRE (NP_000374.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - AIRE Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using AIRE antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - AIRE Polyclonal Antibody

Immunofluorescence analysis of A549 cells using AIRE antibody . Blue: DAPI for nuclear staining.