

Phospho-GRIN1-S890 pAb



Catalog Number: 166532

Product name

Phospho-GRIN1-S890 pAb

Specificity

Human, Rat

Antibody description

Polyclonal antibody to Phospho-GRIN1-S890

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S890 of human GRIN1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-GRIN1-S890 pAb

Western blot analysis of extracts of C6 cells, using Phospho-GRIN1-S890 antibody at 1:1000 dilution. C6 cells were treated by PMA/TPA (200nM) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 90s.