

GEF-H1 antibody



Catalog Number: 110932

Product name

GEF-H1 antibody

Specificity

Human; other species not tested.

Antibody description

GEF-H1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 111 kDa

Preparation

This antibody was obtained by immunization of GEF-H1 recombinant protein (Accession Number: NM_004723). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110932(GEF-H1 Antibody) at dilution of 1:1000