

Anti-CTDSP1 antibody



Catalog Number: 175932

Product name

Anti-CTDSP1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to CTDSP1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ctdsp1. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.; Western Blot analysis of various cells using CTDSP1 Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).