

Anti-hHR23A Rabbit antibody



Catalog Number: 167932

Product name

Anti-hHR23A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

hHR23A Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human hHR23A

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/2000-1/10000

FC: 1/20

Validations



Western blot detection of hHR23A in HeLa, A549, HL-60, U2OS using hHR23A Rabbit mAb(1:1000 diluted)