

# ELOVL6 antibody

Catalog Number: 110232

## Product name

ELOVL6 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ELOVL6 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human hepat°Cirrhosis. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of ELOVL6 recombinant protein (Accession Number: XM\_011532233). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:20-1:200

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:110232(ELOVL6 antibody) at dilution of 1:800 

Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:110232(ELOVL6 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:110232(ELOVL6 antibody) at dilution of 1:100 (under 40x lens)