

Anti-FCRLA antibody



Catalog Number: 141232

Product name

Anti-FCRLA antibody

Specificity

Human FCRLA

Antibody description

Rabbit Polyclonal to Human FCRLA

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human FCRLA, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Items	Lanes	A	B
Sample (whole cell lysate)		Raji	Hela
Sample Volume (µg/lane)		30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	5-20 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.		
Explanation	Developed using Odyssey imaging system. Predicted band size : 39 kDa Observed band size : 48 kDa		

Human FCRLA Western blot (WB) 14874