

GTF2H1 antibody



Catalog Number: 111232

Product name

GTF2H1 antibody

Specificity

human, mouse; other species not tested.

Antibody description

GTF2H1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 62kd

Preparation

This antibody was obtained by immunization of GTF2H1 recombinant protein (Accession Number: NM_005316). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111232(GTF2H1 antibody) at dilution of 1:200