

MEF2C antibody

Catalog Number: 112632

Product name

MEF2C antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

MEF2C Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse colon tissue. Positive IP detected in SH-SY5Y cells. Positive IHC detected in human lymphoma tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 52 kDa; 63kDa

Preparation

This antibody was obtained by immunization of MEF2C recombinant protein (Accession Number: XM_047417182). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:112632(MEF2C antibody) at dilution of 1:500



IP Result of anti-MEF2C (IP:Catalog No:112632, 3ug; Detection:Catalog No:112632 1:500) with SH-SY5Y cells lysate 2000ug.




MEF2C antibody



Catalog Number: 112632

Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:112632(MEF2C antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:112632(MEF2C Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)