

RFTN2 antibody

Catalog Number: 114632

Product name

RFTN2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RFTN2 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse brain tissue, mouse testis tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 56-60 kDa

Preparation

This antibody was obtained by immunization of RFTN2 recombinant protein (Accession Number: NM_144629). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:400


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114632(RFTN2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:114632(RFTN2 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:114632(RFTN2 Antibody) at dilution of 1:100

RFTN2 antibody



Catalog Number: 114632

(under 40x lens).