

Anti-IL13RA1/IL-13RA1 antibody



Catalog Number: 105632

Product name

Anti-IL13RA1/IL-13RA1 antibody

Immunogen

[Human IL13RA1/IL-13RA1 \(His Tag\) recombinant protein](#)

Specificity

Human IL13RA1 / IL-13RA1

Antibody description

Rabbit monoclonal to IL13RA1/IL-13RA1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human IL13RA1 / IL-13RA1 (rh IL13RA1 / IL-13RA1; NP_001551.1; Met1-Thr343) and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of

water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



IL13RA1 / IL-13RA1 Antibody (APC), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of anti-human IL13RA1(CD213a1) on MCF-7 cells. MCF-7 cells were stained with APC-anti human IL13RA1(CD213a1)(solid line) or APC-rabbit isotype control (dashed line).