

TP63 antibody



Catalog Number: 107632

Product name

TP63 antibody

Specificity

Human; other species not tested.

Antibody description

TP63 Mouse Monoclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of TP63 recombinant protein (Accession Number: NM_003722). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107632(TP63 Antibody) at dilution of 1:1000