

Human ATG5/APG5L/ASP recombinant protein



Catalog Number: 519632

General Information

Gene Name Synonym

Autophagy protein 5, APG5-like, Apoptosis-specific protein, APG5L, ASP, ATG5,

Protein Construction

A DNA sequence encoding the Human ATG5/APG5L/ASP Met1-Asp275 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 0.2M NaCl, 1mM DTT, 40% glycerol, pH 8.0 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ATG5/APG5L/ASP recombinant protein