

SLC25A3 antibody



Catalog Number: 115332

Product name

SLC25A3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SLC25A3 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, MCF7 cells. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of SLC25A3 recombinant protein (Accession Number: NM_002635). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:115332(SLC25A3 antibody) at dilution of 1:500