

Anti-hptII (2D3) Mouse antibody



Catalog Number: 168332

Product name

WB

Anti-hptII (2D3) Mouse antibody

Dilutions

Antibody description

WB: 1/1000

hptII (2D3) Mouse monoclonal antibody

Validations

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications



Western blot detection of HptII in Transgenic rice (leaf) and Non-transgenic rice (leaf) cell lysates using HptII mouse mAb (1:1000 diluted). Predicted band size: 36, 38KDa. Observed band size: 36, 38KDa. Beta-actin mouse mAb(200068-8F10, 1:1000 diluted) was used as a loading control.