

PEX7 antibody



Catalog Number: 113732

Product name

PEX7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PEX7 Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of PEX7 recombinant protein (Accession Number: NM_000288). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:113732(PEX7 antibody) at dilution of 1:800