

VDAC3 antibody



Catalog Number: 116732

Product name

VDAC3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

VDAC3 Rabbit Polyclonal antibody. Positive IHC detected in human osteosarcoma. Positive WB detected in human heart tissue, human brain tissue. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of VDAC3 recombinant protein (Accession Number: NM_005662). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000


IHC: 1:20-1:200

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116732(VDAC3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human osteosarcoma using Catalog No:116732(VDAC3 antibody) at dilution of 1:100 (under 10x lens) 

Immunohistochemical of paraffin-embedded human osteosarcoma using Catalog No:116732(VDAC3 antibody) at dilution of 1:100 (under 40x lens)