

CTRP1 antibody

Catalog Number: 109732

Product name

CTRP1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CTRP1 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human placenta tissue. Positive IF detected in HepG2 cells. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of CTRP1 recombinant protein (Accession Number: NM_030968). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109732(C1QTNF1 antibody) at dilution of 1:500




Immunofluorescent analysis of HepG2 cells, using C1QTNF1 antibody Catalog No:109732 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



CTRP1 antibody

Catalog Number: 109732

Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:109732(C1QTNF1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:109732(C1QTNF1 Antibody) at dilution of 1:50 (under 40x lens)