

# NCOA3 antibody

Catalog Number: 113042

## Product name

NCOA3 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


## Specificity

Human; other species not tested.

## Validations

## Antibody description

NCOA3 Rabbit Polyclonal antibody. Positive FC detected in Raji cells. Positive IF detected in HepG2 cells. Positive IHC detected in human lung cancer tissue, human heart tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 160kd



## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006534). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113042(NCOA3 antibody) at dilution of 1:300

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IF, IHC, FC


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunofluorescent analysis of HepG2 cells, using NCOA3 antibody Catalog No:113042 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).




Immunohistochemistry of paraffin-embedded


# NCOA3 antibody

Catalog Number: 113042

human lung cancer tissue slide using Catalog No:113042(NCOA3 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:113042(NCOA3 Antibody) at dilution of 1:50



1X10<sup>6</sup> Raji cells were stained with 0.2ug NCOA3 antibody (Catalog No:113042, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.