

Anti-KDEL2 antibody



Catalog Number: 176042

Product name

Anti-KDEL2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to KDEL2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human erd22. aa range:81-130

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations



Fig1: Westerns Blot analysis of various cells using KDEL Receptor 2 Polyclonal Antibody diluted at 1:1000