

ZFP36L2 antibody



Catalog Number: 117042

Product name

ZFP36L2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZFP36L2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HL-60 cells. Observed molecular weight by Western-blot: 51kd; 60kd

Preparation

This antibody was obtained by immunization of ZFP36L2 recombinant protein (Accession Number: NM_006887). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117042(ZFP36L2 antibody) at dilution of 1:500