

Anti-DAZAP1 Rabbit antibody



Catalog Number: 169042

Product name

Anti-DAZAP1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

DAZAP1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human DAZAP1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

1:2000-1:10000

1:20

1:20

Validations



Western blot detection of DAZAP1 in K562,C6,3T3,Hela cell lysates using DAZAP1 Rabbit mAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.