

HOXD8 antibody



Catalog Number: 111442

Product name

HOXD8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HOXD8 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of HOXD8 recombinant protein (Accession Number: NM_001199746). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:111442(HOXD8 Antibody) at dilution of 1:300