

P21 antibody

Catalog Number: 107442

Product name

Monoclonal

P21 antibody

Ig Type

Specificity

Mouse IgG2b

Human, Mouse, Monkey; other species not tested.

Applications

Antibody description

ELISA, WB

P21 Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells. Observed molecular weight by Western-blot: 21 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of P21 recombinant protein (Accession Number: NM_001220777). Purification method: Protein A purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107442(P21 Antibody) at dilution of 1:1000