

Anti-beta 2 Microglobulin (4C3) Mouse antibody



Catalog Number: 168442

Product name

Anti-beta 2 Microglobulin (4C3) Mouse antibody

Specificity

Human

Antibody description

beta 2 Microglobulin (4C3) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human β 2-microglobulin protein expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b

Applications

ELISA

Validations



sELISA standard Curve for u03b22-microglobulin.