

# GALK1 antibody



Catalog Number: 110842

---

## Product name

GALK1 antibody

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

GALK1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, mouse heart tissue. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of GALK1 recombinant protein (Accession Number: NM\_001381985). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110842(GALK1 antibody) at dilution of 1:1200