

GPRIN3-Specific antibody



Catalog Number: 111142

Product name

GPRIN3-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GPRIN3-Specific Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 82 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_198281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111142(GPRIN3-Specific antibody) at dilution



of 1:300

Immunofluorescent analysis of HepG2 cells, using GPRIN3 antibody Catalog No:111142 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).