

# PPP1CB antibody

Catalog Number: 114142

## Product name

PPP1CB antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPP1CB Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in A549 cells, mouse brain tissue. Observed molecular weight by Western-blot: 37-40 kDa or 33 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002709). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:50-1:500

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:114142(PPP1CB antibody) at dilution of 1:2000



Immunofluorescent analysis of (10% Formaldehyde) fixed A549 cells using Catalog No:114142(PPP1CB Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)