

Ozuriftamab reference antibody



Catalog Number: 177142

Product name

Ozuriftamab reference antibody

Specificity

Human ROR2

Antibody description

Anti-ROR2 Antibody (Ozuriftamab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-ROR2 Antibody (Ozuriftamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-ROR2 Antibody (Ozuriftamab) is 99.73%, determined by SEC-HPLC.



Bioactivity: ELISA

Immobilized human ROR2 His at 2 ug/mL can bind Anti-ROR2 Antibody (Ozuriftamab), EC50=0.114 ug/mL.