

# HSD11B1 antibody



Catalog Number: 111542

---

## Product name

HSD11B1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HSD11B1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of HSD11B1 recombinant protein (Accession Number: NM\_005525). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111542(HSD11B1 antibody) at dilution of 1:400